

FIG. 1

009780" 0696E960




FIG. 2

009T80" 0696E960

100

FOOD SUPPLY, PROCESS LINE

SAMPLE

PREPARE CULTURE




FOR PCR

mp 30

DATA MINING SOFTWARE

DATA BASES

DISPLAY



PROCESSED ARRAY

ROBOT, 30

LOAD, PROCESS




ARRAY, 50

PROCESS CONTROL KNOWLEDGE/ FEEDBACK

NEW ARRAY KNOWLEDGE/ FEEDBACK

FIG. 3



Gene target table for model array

	A-B	C-D	E-F	G-H	I-J	K-L	M-N	O-P	Q-R	S-T	U-V	W-X
1	Bacillus Hemolysin BL (hblA)	Bacillus Enterotoxin B (entB)	Bacillus dacF	Bacillus spoIIIGA	Bacillus sigE	Bacillus sigF	Bacillus Enterotoxin A (hblB)	Bacillus spoOA	Bacillus spoIIAA	Listeria ipa-p60	Listeria ltaA	Listeria mpl
2	Blank	Blank	Control gfp	Control gfp 1:10	Control gfp 1:100	Control gfp 1:1000			Control gfp	Control gfp		
3	Campylobacter/Listeria	Campylobacter cipA	Campylobacter lig3	Campylobacter flaB	Campylobacter ciaE	Campylobacter folP	Campylobacter pglD	Campylobacter pglF	Campylobacter pglB	Listeria InlB	Listeria hlyA	Listeria InlA
4	E.coli/Listeria	E.coli aspRNA	E.coli lacZ	E.coli thrRNA	E.coli hlyD	E.coli gdh	E.coli rpoS	E.coli uidA	E.coli air	Listeria prfA	Listeria plcB	Listeria plcA
5	Yersinia/O157	Yersinia yopB	Yersinia yopH	Yersinia yopE	Yersinia lcrV	Yersinia virF	Yersinia yopN	Yersinia yopP	Yersinia yopM	O157:H7 ehxA	O157:H7 hlyA	O157:H7 toxB
6	Staphylococcus aureus/E. coli O157:H7	S. aureus sarA	S. aureus entA	S. aureus entD	S. aureus nuclease	S. aureus set1A	S. aureus agrB	S. aureus entI	S. aureus tsst-1	O157:H7 espP	O157:H7 eeeA	O157:H7 elpC
7	Shigella/E. coli O157:H7	Shigella gfp	Shigella lcsB	Shigella stxA	Shigella ipaC	Shigella vacB	Control gfp	Shigella ipaB	Shigella virB	Shigella sod	Control gfp	Shigella virF
8	Salmonella/E. coli O157:H7	Salmonella pagC	Salmonella invJ	Salmonella spiA	Salmonella prgK	Salmonella traJ	Salmonella traK	Salmonella invH	Salmonella traB	O157:H7 Rof2	O157:H7 stx2	O157:H7 escT

FIG. 4

009T80" 0696E960

1A-1R = *Bacillus cereus*
2A-2Z = controls
3A-3R = *Campylobacter jejuni*
4A-4R = *E. coli*
5A-5R = *Yersinia enterocolitica*
6A-6R = *Staphylococcus aureus*
7A-B, 7E-L, 7O-T, 7W-X = *Shigella dysenteriae*
8A-8R = *Salmonella enterica*

1S-1X = *Listeria monocytogenes*
2S-2X = controls (UV are positive controls)
3S-3X = *Listeria monocytogenes*
4S-4X = *Listeria monocytogenes*
5S-5X = *E. coli* O157:H7
6S-6X = *E. coli* O157:H7
7C-D, 7M-N, 7U-V = controls (UV are + controls)
8S-8X = *E. coli* O157:H7

