

FIG. 1

009780" 0696E960




FIG. 2

100

FOOD SUPPLY, PROCESS LINE

SAMPLE

PREPARE CULTURE

SOLVENTS
BUFFERS
PRIMERS, ET

DATA MINING
SOFTWARE

mp 30

DATA BASES

DISPLAY

READER
60

PROCESSED
ARRAY

ROBOT, 30


LOAD, PROCESS

ARRAY, 50

PROCESS CONTROL
KNOWLEDGE/
FEEDBACK

NEW ARRAY
KNOWLEDGE/
FEEDBACK

FIG. 3



Gene target table for model array

	A-B	C-D	E-F	G-H	I-J	K-L	M-N	O-P	Q-R	S-T	U-V	W-X
1	Bacillus Hemolysin BL (hblA)	Bacillus Enterotoxin B (entB)	Bacillus dacF	Bacillus spoIIIGA	Bacillus sigE	Bacillus sigF	Bacillus Enterotoxin A (hblB)	Bacillus spoOA	Bacillus spoIIAA	Listeria ipa-p60	Listeria ltaA	Listeria mpl
2	Blank	Blank	Control gfp	Control gfp 1:10	Control gfp 1:100	Control gfp 1:1000			Control gfp	Control gfp		
3	Campylobacter/Listeria	Campylobacter cipA	Campylobacter lig3	Campylobacter flaB	Campylobacter ciaE	Campylobacter folP	Campylobacter pglD	Campylobacter pglF	Campylobacter pglB	Listeria InlB	Listeria hlyA	Listeria InlA
4	E.coli/Listeria	E.coli asprRNA	E.coli lacZ	E.coli thrRNA	E.coli hlyD	E.coli gdh	E.coli rpoS	E.coli uidA	E.coli alr	Listeria prfA	Listeria plcB	Listeria plcA
5	Yersinia/O157	Yersinia yopD	Yersinia yopH	Yersinia yopE	Yersinia lcrV	Yersinia virF	Yersinia yopN	Yersinia yopP	Yersinia yopM	O157:H7 ehxA	O157:H7 hlyA	O157:H7 toxB
6	Staphylococcus aureus/E. coli O157:H7	S. aureus sarA	S. aureus entA	S. aureus entD	S. aureus nuclease	S. aureus sar1A	S. aureus agrB	S. aureus entI	S. aureus tsst-1	O157:H7 espP	O157:H7 eeeA	O157:H7 elpC
7	Shigella/E. coli O157:H7	Shigella fur	Shigella lcsB	Shigella stxA	Shigella ipaC	Shigella vacB	Control gfp	Shigella ipaB	Shigella virB	Shigella sod	Control gfp	Shigella virF
8	Salmonella/E. coli O157:H7	Salmonella pagC	Salmonella invJ	Salmonella spiA	Salmonella prgK	Salmonella traJ	Salmonella traK	Salmonella invH	Salmonella traB	O157:H7 Rof2	O157:H7 stx2	O157:H7 escT

FIG. 4

009T80" 0696E960

1A-1R = *Bacillus cereus*
 2A-2Z = controls
 3A-3R = *Campylobacter jejuni*
 4A-4R = *E. coli*
 5A-5R = *Yersinia enterocolitica*
 6A-6R = *Staphylococcus aureus*
 7A-B, 7E-L, 7O-T, 7W-X = *Shigella dysenteriae*
 8A-8R = *Salmonella enterica*

1S-1X = *Listeria monocytogenes*
 2S-2X = controls (UV are positive controls)
 3S-3X = *Listeria monocytogenes*
 4S-4X = *Listeria monocytogenes*
 5S-5X = *E. coli* O157:H7
 6S-6X = *E. coli* O157:H7
 7C-D, 7M-N, 7U-V = controls (UV are + controls)
 8S-8X = *E. coli* O157:H7





FIG 4A